

ORP (redox) measurement. Metallic electrodes. A little theory

ORP Redox potential

Oxide-reduction potential, ORP, also known as redox potential, indicates the oxidising or reducing properties of a substance.

To measure ORP, a pH-meter which measures mV and a metallic electrode made of platinum or gold are used.

Platinum is recommended for highly-oxidising solutions that contain chlorides, redox titrations, etc.

Gold is recommended for highly acidic solutions,

titrations of ampicillin or penicillin with mercury nitrate, for redox titrations where chrome or iron etc, are present.

Silver electrodes are also available.

Silver is mainly used in argentometric titration.

Calibration of metallic electrodes

Metallic electrodes do not present significant drift of potential, which is why they are not calibrated.

However, some deviations of potential may occur after continued use of the electrode as a result of alterations of the metallic surface or contamination of the reference electrode. Redox standard solutions are available to check that these electrodes are working well, see page 53.

Effects of temperature on ORP (redox) measuring...

... on the electrode

Temperature has no significant effect whatsoever on metal electrodes.

... on a given sample

The oxidation-reduction potential of a sample depends on temperature. However, pH-meters do not provide any type of compensation.

In this case, connecting a temperature sensor to the pH-meter will serve only to find out the temperature of the sample.

Combined metallic electrodes Essential components

Insofar as their connectors, reference elements, electrolyte and diaphragms are concerned, combined metallic electrodes are the same as pH electrodes. However, the glass membrane is replaced with a ring or wire made from a noble metal.



Connector

With S7 screw cap, to attach the cable.

Body material

Normally made of glass.

Reference element

This is a cell that supplies a stable potential (mV). There are a number of types:

Silver wire (Ag)

Galvanically coated with AgCl, it is the typical reference element.

Encapsulated AgCl crystals

The silver wire comes into contact with a portion of the AgCl crystals inside a small glass tube.

Electrolyte

This is a concentrated saline solution in contact with the reference element. Metallic electrodes are available with liquid electrolytes.

Diaphragm

This is the point of union between the electrolyte and the sample. It is the critical part of the electrode as it has a direct influence on its life-span.

A wide range of diaphragms are available. They vary according to manufacturer, application, quality of electrode and price.

This catalogue features electrodes with sleeve, ceramic and open diaphragms.



The sleeve diaphragm consists of a hole in the body of the electrode, half-closed by a PTFE ring which provides a high level of electrolyte flow..



The ceramic diaphragm is a chemically-inert, porous ceramic plate. These are the traditional diaphragms. It allows a tiny amount of electrolyte to flow through.

Metal wire or ring

May be made of either silver, platinum or gold.

Practical considerations

Immersion depth

To obtain a correct measurement of the ORP, the electrode must be immersed until the diaphragm is covered.

Minimum volume of sample

This will depend on the shape of the measuring recipient.

Price-quality ratio

In practice, the quality of the electrode used has a direct bearing on the reliability of a measurement. CRISON only supplies top quality electrodes because experience has proven them to be a better investment in the medium and long-term.

Common problems

Use of the wrong electrode, or mis-use of the electrode, causes a series of problems that will substantially reduce the lifespan of the electrode. The most common problems, along with their causes and consequences, are:

- Diaphragm blockage.
Open measuring circuit. Unstable readings...
- Reference system contamination.
Drift of mV values.
- Build-up of dirt on metallic surface.
Slow response and wrong measurements.
- Slow response.
As described above.

Guarantee

CRISON electrodes are guaranteed for a six-month period.

The guarantee only covers manufacturing defects and does not cover defects that could arise from incorrect use, handling, application or maintenance, or as a result of premature wear inherent to certain samples.

ORP

Applications for metallic electrodes

Application	Electrode	Metal	Electrolyte	Observations
General ORP measurement	52 61	Platinum	CRISOLYT A	The 50 55 model, especially for the PH 25. This type of sample may contain elements which contaminate the diaphragm, which is why CRISON recommends electrodes with a clog-resistant diaphragm.
...in waste water	50 55	Platinum	CRISOLYT	
... in galvanizing baths	52 62 / 50 55	Platinum	CRISOLYT	
Redox titrations in general	52 63	Gold	CRISOLYT	These are special-length electrodes, which makes them appropriate for taking measurements directly in COD tubes. Mainly used with CRISON'S special COD samplers (Sampler 20).
...COD analysis	52 66	Platinum	CRISOLYT A	
...in microsamples	52 70	Gold	CRISOLYT A	
...in the presence of chrome or iron	52 65	Platinum	CRISOLYT A	
Argentometric titrations	52 63	Gold	CRISOLYT A	Gold is recommended as an indicator.
Argentometry in microsamples	52 69	Gold	CRISOLYT	
Argentometric titrations	52 60	Silver	KNO ₃ 1M	An electrode with silver as an indicator metal should be used. In the case of the micro electrode, it can be used to work on microanalysis with the CRISON Sampler 45.
Argentometry in microsamples	52 68	Silver	KNO ₃ 1M	
KF Titrations	52 64	Double platinum	-	The technique used in this type of analysis is bipotentiometry. They are indicator electrodes. The presence of a reference electrode is not required.
SO ₂ Analysis	52 64	Double platinum	-	