




BALANCES SERIES

4Y SERIES



Display

5.7" resistive colour touch screen

Keyboard

On-screen + 8 keys

Customized Display

Workspace + Toolbar

Databases

19

Proximity Sensors

2 programmable IR sensors

Level system

Electronic

Interfaces

2×USB-A, 2×RS 232, Ethernet,
4×IN, 4×OUT, Wireless Connection

X2 SERIES



Display

5" capacitive colour touch screen

Keyboard

On-screen + 6 keys

Customized Display

Workspace + Toolbar

Databases

7

Proximity Sensors

2 programmable IR sensors

Level system

Bull's eye level

Interfaces

USB-A, USB-B, 2×RS 232,
Ethernet, Wireless Connection

R SERIES



Display

LCD backlit

Keyboard

Membrane, 14 keys

Customized Display

–

Databases

5

Proximity Sensors

–

Level system

Bull's eye level

Interfaces

USB-A, USB-B, 2×RS 232,
Ethernet, Wireless Connection