

HOM-302SS




**HOM-302SS/312SS, Micronizing Homo Mixer**

**Applications:**

This machine uses the pressurized, accelerated, destroyed andstricken interaction between the liquid molecules flowing in, the fluid particles are even to spout out instantly by the side of the gaps, and then through multiple shearing, grinding, dispersing to damage into fine and integration and micronized effect. Offering a great help to improve the quality of foodstuffs (like dairy products, fruitjuices, vinegar, oil and syrup) or paint, chemicals and various products; the model can emulsifying the products, and not precipitated, not isolated, and getting meticulous are the biggest features.

**Features:**

- Spring-lift pole for easy operation
- HOM-302SS is with digital speed display
- HOM-312SS is with time & speed displays, and digital time setting range 0-99H 0-59Min.
- Impeller is made of stainless steel SUS316
- Stepless & stable speed; auto torque compensated function
- With the over heat & fuse protection devices.



Model	HOM-302SS	HOM-312SS
Motor	DC 220W	
Stirrer Speed (rpm)	500 ~ 12,000	
Capacity (L)	2 (H <sub>2</sub> O)	
Timer	-	0 ~ 99 Hr, 59 Min
Weight (kg)	10	
Body Size (W x L x H mm)	166 x 283 x 415	
Overall dimension (W x L x H mm)	328 x 303 x 540	
Power	AC 110V/220V, 50/60Hz	